




13485



# Spinal Solutions



[www.4umedical.com.tr](http://www.4umedical.com.tr)



**Medical**

A full-page background image showing a person in a dynamic, athletic pose climbing a dark, craggy rock face. The person is illuminated from below, creating a strong silhouette effect against the dark rock. A large, bright, circular light source, possibly the moon, is visible in the upper right corner of the sky.

Keeps You  
Moving



**PLEASE SCAN QR CODE  
TO LEARN MORE!**



# Contents

• Posterior Cervical Spinal Screw System

• Anterior Cervical Plate

• Cervical Corpectomy Cage

• Cervical Corpectomy Plate Cage

• Cervical Disc Prosthesis

• Cervical PEEK Cages

• Posterior Spinal Screw System

• Posterior Spinal Rods

• Thoracolumbar Connectors

• Thoracolumbar Hooks

• Thoracolumbar Corpectomy Cage

• Thoracolumbar Interspinous Device

• Thoracolumbar PLIF & TLIF PEEK Cages

• Instrument Sets



# Introduction



The 4U Spinal Solutions offers a comprehensive portfolio for neurosurgeons specializing in spinal fusion. The system is designed to provide procedural solutions for the cervical, thoracic, lumbar and/or sacral spinal segments.

The 4U System is used to treat defects and stabilize the spine, joint or other skeletal structures. 4U Spinal Fixation Systems are the next generation of top-loading, top-tightening systems designed to spine for optimal deformity and correction. This requires a comprehensive selection of implants and well-designed instrumentation.

Therefore, 4U Spinal Solutions offers a wide range and high quality high quality range of implants that support anterior & posterior correction techniques. It consists of long term surgical and implantable devices, different sizes of polyaxial self-tapping pedicle screws in various designs, rods, connectors, hooks and different types of interbody cages.

# Posterior Cervical Spinal Screw System



## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- Angled & Anatomical design
- Easier screw fixation
- Flexibility for screw and rod replacement

### 4U SPINAL POSTERIOR CERVICAL POLYAXIAL SCREW

| Code    | Size      | Code    | Size      |
|---------|-----------|---------|-----------|
| 1803010 | 3.0x10 mm | 1804010 | 4.0x10 mm |
| 1803012 | 3.0x12 mm | 1804012 | 4.0x12 mm |
| 1803014 | 3.0x14 mm | 1804014 | 4.0x14 mm |
| 1803016 | 3.0x16 mm | 1804016 | 4.0x16 mm |
| 1803018 | 3.0x18 mm | 1804018 | 4.0x18 mm |
| 1803020 | 3.0x20 mm | 1804020 | 4.0x20 mm |
| 1803022 | 3.0x22 mm | 1804022 | 4.0x22 mm |
| 1803024 | 3.0x24 mm | 1804024 | 4.0x24 mm |
| 1803026 | 3.0x26 mm | 1804026 | 4.0x26 mm |
| 1803028 | 3.0x28 mm | 1804028 | 4.0x28 mm |
| 1803030 | 3.0x30 mm | 1804030 | 4.0x30 mm |
| 1803032 | 3.0x32 mm | 1804032 | 4.0x32 mm |
| 1803510 | 3.5x10 mm | 1804510 | 4.5x10 mm |
| 1803512 | 3.5x12 mm | 1804512 | 4.5x12 mm |
| 1803514 | 3.5x14 mm | 1804514 | 4.5x14 mm |
| 1803516 | 3.5x16 mm | 1804516 | 4.5x16 mm |
| 1803518 | 3.5x18 mm | 1804518 | 4.5x18 mm |
| 1803520 | 3.5x20 mm | 1804520 | 4.5x20 mm |
| 1803522 | 3.5x22 mm | 1804522 | 4.5x22 mm |
| 1803524 | 3.5x24 mm | 1804524 | 4.5x24 mm |
| 1803526 | 3.5x26 mm | 1804526 | 4.5x26 mm |
| 1803528 | 3.5x28 mm | 1804528 | 4.5x28 mm |
| 1803530 | 3.5x30 mm | 1804530 | 4.5x30 mm |
| 1803532 | 3.5x32 mm | 1804532 | 4.5x32 mm |

## Posterior Cervical Spinal System



### Cervical Rod 3.5


| Code    | Size       |
|---------|------------|
| 2400330 | 3.5x30 mm  |
| 2400340 | 3.5x40 mm  |
| 2400350 | 3.5x50 mm  |
| 2400360 | 3.5x60 mm  |
| 2400370 | 3.5x70 mm  |
| 2400380 | 3.5x80 mm  |
| 2400390 | 3.5x90 mm  |
| 2403100 | 3.5x100 mm |
| 2403120 | 3.5x120 mm |
| 2403150 | 3.5x150 mm |
| 2403180 | 3.5x180 mm |
| 2403200 | 3.5x200 mm |
| 2403240 | 3.5x240 mm |

### Cervical Lateral Connector

| Code    | Size  |
|---------|-------|
| 2700010 | 10 mm |
| 2700018 | 18 mm |



# Posterior Cervical Spinal Screw System



Posterior Cervical  
Spinal System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- Angled & Anatomical design
- Easier screw fixation
- Flexibility for screw and rod replacement

### Multiaxial Transverse Connector

| Code    | Size  |
|---------|-------|
| 2600030 | 30 mm |
| 2600035 | 35 mm |
| 2600040 | 40 mm |
| 2600045 | 45 mm |
| 2600050 | 50 mm |
| 2600055 | 55 mm |

### Occipital Plate

| Code    | Size |
|---------|------|
| 2450005 | L    |
| 2450010 | M    |
| 2450015 | S    |

### Cervical Hook

| Code    | Size |
|---------|------|
| 3000044 | 4X4  |
| 3000045 | 4X5  |
| 3000046 | 4X6  |
| 3000047 | 4X7  |



### Transverse Connector

| Code    | Size  |
|---------|-------|
| 2500030 | 30 mm |
| 2500035 | 35 mm |
| 2500040 | 40 mm |
| 2500045 | 45 mm |



### Occipital Plate Screw

| Code    | Size     |
|---------|----------|
| 2454004 | 4.0x4 mm |
| 2454005 | 4.0x5 mm |
| 2454006 | 4.0x6 mm |
| 2454504 | 4.5x4 mm |
| 2454505 | 4.5x5 mm |
| 2454506 | 4.5x6 mm |



### Domino-Connector

#### Code

27000050  
27000060




# Anterior Cervical Plate

## Features & Benefits

- Ti6Al4V-ELI raw material
- Easy locking system tighten screws
- Low profile
- Anatomical design





Length(L)



### Cervical Plate

### Cervical Plate Screws

| Code    | Size  | Level | Code    | Size      |
|---------|-------|-------|---------|-----------|
| 8000017 | 17 mm | 1     | 8103512 | 3.5x12 mm |
| 8000020 | 20 mm | 1     | 8103514 | 3.5x14 mm |
| 8000023 | 23 mm | 1     | 8103516 | 3.5x16 mm |
| 8000025 | 25 mm | 1     | 8103518 | 3.5x18 mm |
| 8000027 | 27 mm | 2     | 8103520 | 3.5x20 mm |
| 8000030 | 30 mm | 2     |         |           |
| 8000033 | 33 mm | 2     | 8104012 | 4.0x12 mm |
| 8000036 | 36 mm | 2     | 8104014 | 4.0x14 mm |
| 8000040 | 40 mm | 2     | 8104016 | 4.0x16 mm |
| 8000045 | 45 mm | 3     | 8104018 | 4.0x18 mm |
| 8000050 | 50 mm | 3     | 8104020 | 4.0x20 mm |
| 8000055 | 55 mm | 3     |         |           |
| 8000060 | 60 mm | 3     | 8104512 | 4.5x12 mm |
| 8000065 | 65 mm | 3     | 8104514 | 4.5x14 mm |
| 8000070 | 70 mm | 4     | 8104516 | 4.5x16 mm |
| 8000075 | 75 mm | 4     | 8104518 | 4.5x18 mm |
| 8000080 | 80 mm | 4     | 8104520 | 4.5x20 mm |
| 8000085 | 85 mm | 4     |         |           |
| 8000090 | 90 mm | 4     |         |           |




Anterior Cervical Plate

# Cervical Corpectomy Cage

## Features & Benefits

- Ti6Al4V-ELI raw material
- One-easy locking mechanism and teethed and angled superior and inferior surface minimize the displacement and coming out risk
- Rotatable superior endplate
- Efficient graft space
- Anatomical design
- Angled superior surface option



| Code    | D  | CL | OL |
|---------|----|----|----|
| 6001013 | 10 | 10 | 13 |
| 6001017 | 10 | 13 | 17 |
| 6001625 | 10 | 16 | 25 |
| 6101013 | 12 | 10 | 13 |
| 6101017 | 12 | 13 | 17 |
| 6101625 | 12 | 16 | 25 |
| 6102440 | 12 | 24 | 40 |
| 6103965 | 12 | 39 | 65 |
| 6201013 | 14 | 10 | 13 |
| 6201017 | 14 | 13 | 17 |
| 6201625 | 14 | 16 | 25 |
| 6202440 | 14 | 24 | 40 |
| 6203965 | 14 | 39 | 65 |
| 6301013 | 16 | 10 | 13 |
| 6301017 | 16 | 13 | 17 |
| 6301625 | 16 | 16 | 25 |
| 6302440 | 16 | 24 | 40 |
| 6303965 | 16 | 39 | 65 |
| 6401013 | 20 | 10 | 13 |
| 6401017 | 20 | 13 | 17 |
| 6401625 | 20 | 16 | 25 |
|         |    |    |    |
| 6151625 | 12 | 16 | 25 |
| 6152440 | 12 | 24 | 40 |
| 6153965 | 12 | 39 | 65 |
| 6251625 | 14 | 16 | 25 |
| 6252440 | 14 | 24 | 40 |
| 6253965 | 14 | 39 | 65 |




## Angled

# Cervical Corpectomy Plate Cage



## Features & Benefits

- Ti6Al4V-ELI raw material
- One-easy locking mechanism and teethed and angled superior and inferior surface minimize the displacement and coming out risk
- Rotatable superior endplate
- Efficient graft space
- Anatomical design
- Angled superior surface option



Cervical Corpectomy Plate Cage

| Code    | D  | CL | OL |
|---------|----|----|----|
| 6011013 | 10 | 10 | 13 |
| 6011017 | 10 | 13 | 17 |
| 6011625 | 10 | 16 | 25 |
| 6111013 | 12 | 10 | 13 |
| 6111017 | 12 | 13 | 17 |
| 6111625 | 12 | 16 | 25 |
| 6112440 | 12 | 24 | 40 |
| 6113965 | 12 | 39 | 65 |
| 6211013 | 14 | 10 | 13 |
| 6211017 | 14 | 13 | 17 |
| 6211625 | 14 | 16 | 25 |
| 6212440 | 14 | 24 | 40 |
| 6213965 | 14 | 39 | 65 |
| 6311013 | 16 | 10 | 13 |
| 6311017 | 16 | 13 | 17 |
| 6311625 | 16 | 16 | 25 |
| 6312440 | 16 | 24 | 40 |
| 6313965 | 16 | 39 | 65 |
| 6411013 | 20 | 10 | 13 |
| 6411017 | 20 | 13 | 17 |
| 6411625 | 20 | 16 | 25 |



# Cervical Disc Prosthesis




## Cervical Disc Prosthesis

### Features & Benefits

- Ti6Al4V-ELI raw material
- Toothed superior & inferior surfaces offers appropriate fixation between endplates
- ( $\pm$ ) 10° extension / flexion and lateral bending, and ( $\pm$ ) 8° rotation of the discs
- Compatible with MRI and CT

| Code    | L  | W  | H |
|---------|----|----|---|
| 7001204 | 12 | 14 | 4 |
| 7001264 | 12 | 16 | 4 |
| 7101404 | 14 | 14 | 4 |
| 7101464 | 14 | 16 | 4 |
| 7001205 | 12 | 14 | 5 |
| 7001265 | 12 | 16 | 5 |
| 7101405 | 14 | 14 | 5 |
| 7101465 | 14 | 16 | 5 |
| 7001206 | 12 | 14 | 6 |
| 7001266 | 12 | 16 | 6 |
| 7101406 | 14 | 14 | 6 |
| 7101466 | 14 | 16 | 6 |
| 7001207 | 12 | 14 | 7 |
| 7001267 | 12 | 16 | 7 |
| 7101407 | 14 | 14 | 7 |
| 7101467 | 14 | 16 | 7 |
| 7001208 | 12 | 14 | 8 |
| 7001268 | 12 | 16 | 8 |
| 7101408 | 14 | 14 | 8 |
| 7101468 | 14 | 16 | 8 |

# Cervical Expandable PEEK Cage



| Code    | L  | W  | H |
|---------|----|----|---|
| 5601205 | 12 | 14 | 5 |
| 5601206 | 12 | 14 | 6 |
| 5601207 | 12 | 14 | 7 |
| 5601208 | 12 | 14 | 8 |


### Features & Benefits

- Poly Ether Ether Ketone (PEEK) and Ti6Al4V-ELI raw material
- Expanding screw offers appropriate expansion
- Titanium Pins Markers and expanding screw helps to be seen in imaging machines
- Compatible with MRI and CT

## Cervical Expandable PEEK Cage

# Cervical PEEK Cage

| Code    | L  | W  | H |
|---------|----|----|---|
| 5400414 | 12 | 14 | 4 |
| 5400416 | 12 | 16 | 4 |
| 5500414 | 14 | 14 | 4 |
| 5500416 | 14 | 16 | 4 |
| 5400514 | 12 | 14 | 5 |
| 5400516 | 12 | 16 | 5 |
| 5500514 | 14 | 14 | 5 |
| 5500516 | 14 | 16 | 5 |
| 5400614 | 12 | 14 | 6 |
| 5400616 | 12 | 16 | 6 |
| 5500614 | 14 | 14 | 6 |
| 5500616 | 14 | 16 | 6 |
| 5400714 | 12 | 14 | 7 |
| 5400716 | 12 | 16 | 7 |
| 5500714 | 14 | 14 | 7 |
| 5500716 | 14 | 16 | 7 |
| 5400814 | 12 | 14 | 8 |
| 5400816 | 12 | 16 | 8 |
| 5500814 | 14 | 14 | 8 |
| 5500816 | 14 | 16 | 8 |




**Cervical PEEK Cage**

## Features & Benefits

- Poly Ether Ether Ketone (PEEK)
- Titanium Pins and toothed superior & inferior surface offer good and anatomical fixation
- Tantalum Markers helps to be seen in Imaging Machines
- Compatible with MRI and CT

# Cervical PEEK Cage with Blade



**Cervical Cage with Blade**


## Features & Benefits

- Poly Ether Ether Ketone (PEEK) and Ti6Al4V-ELI raw material
- Blade offers appropriate fixation between endplates
- Tantalum Markers helps to be seen in imaging machines
- Compatible with MRI and CT

| Code    | L  | W  | H |
|---------|----|----|---|
| 5100414 | 12 | 14 | 4 |
| 5300416 | 12 | 16 | 4 |
| 5200414 | 14 | 14 | 4 |
| 5000416 | 14 | 16 | 4 |
| 5100514 | 12 | 14 | 5 |
| 5300516 | 12 | 16 | 5 |
| 5200514 | 14 | 14 | 5 |
| 5000516 | 14 | 16 | 5 |
| 5100614 | 12 | 14 | 6 |
| 5300616 | 12 | 16 | 6 |
| 5200614 | 14 | 14 | 6 |
| 5000616 | 14 | 16 | 6 |
| 5100714 | 12 | 14 | 7 |
| 5300716 | 12 | 16 | 7 |
| 5200714 | 14 | 14 | 7 |
| 5000716 | 14 | 16 | 7 |
| 5100814 | 12 | 14 | 8 |
| 5300816 | 12 | 16 | 8 |
| 5200814 | 14 | 14 | 8 |
| 5000816 | 14 | 16 | 8 |

# 4U Spinal Cervical Anatomic PEEK Cage Screwed


| Code    | L  | W  | H |
|---------|----|----|---|
| 6400414 | 12 | 14 | 4 |
| 6400416 | 12 | 16 | 4 |
| 6500414 | 14 | 14 | 4 |
| 6500416 | 14 | 16 | 4 |
| 6400514 | 12 | 14 | 5 |
| 6400516 | 12 | 16 | 5 |
| 6500514 | 14 | 14 | 5 |
| 6500516 | 14 | 16 | 5 |
| 6400614 | 12 | 14 | 6 |
| 6400616 | 12 | 16 | 6 |
| 6500614 | 14 | 14 | 6 |
| 6500616 | 14 | 16 | 6 |
| 6400714 | 12 | 14 | 7 |
| 6400716 | 12 | 16 | 7 |
| 6500714 | 14 | 14 | 7 |
| 6500716 | 14 | 16 | 7 |
| 6400814 | 12 | 14 | 8 |
| 6400816 | 12 | 16 | 8 |
| 6500814 | 14 | 14 | 8 |
| 6500816 | 14 | 16 | 8 |



## 4U Spinal Cervical Anatomic PEEK Cage Screwed

### Features & Benefits

- Poly Ether Ether Ketone (PEEK) and Ti6Al4V-ELI raw material
- Toothed superior & inferior surface offer good and anatomical fixation
- Ti6Al4V Titanium Alloy Screws assembled replaces the two spikes(pin) of regular model and anterior cervical plate is not required for supplemental fixation.
- Compatible with MRI and CT



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:
  1. Cannulated screw allowing cement injection
  2. Double-threaded Screw providing a tighter grip
  3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



## Posterior Spinal Screw System

### 4U SPINAL PEDICLE POSTERIOR POLYAXIAL SCREW

| Code    | Size      | Code    | Size      | Code    | Size       | Code    | Size       |
|---------|-----------|---------|-----------|---------|------------|---------|------------|
| 1103520 | 3.5x20 mm | 1105520 | 5.5x20 mm | 1107525 | 7.5x25 mm  | 1108525 | 8.5x25 mm  |
| 1103525 | 3.5x25 mm | 1105525 | 5.5x25 mm | 1107530 | 7.5x30 mm  | 1108530 | 8.5x30 mm  |
| 1103530 | 3.5x30 mm | 1105530 | 5.5x30 mm | 1107535 | 7.5x35 mm  | 1108535 | 8.5x35 mm  |
| 1103535 | 3.5x35 mm | 1105535 | 5.5x35 mm | 1107540 | 7.5x40 mm  | 1108540 | 8.5x40 mm  |
| 1103540 | 3.5x40 mm | 1105540 | 5.5x40 mm | 1107545 | 7.5x45 mm  | 1108545 | 8.5x45 mm  |
| 1103545 | 3.5x45 mm | 1105545 | 5.5x45 mm | 1107550 | 7.5x50 mm  | 1108550 | 8.5x50 mm  |
| 1103550 | 3.5x50 mm | 1105550 | 5.5x50 mm | 1107555 | 7.5x55 mm  | 1108555 | 8.5x55 mm  |
| 1103555 | 3.5x55 mm | 1105555 | 5.5x55 mm | 1107560 | 7.5x60 mm  | 1108560 | 8.5x60 mm  |
| 1104020 | 4.0x20 mm | 1105560 | 5.5x60 mm | 1107570 | 7.5x70 mm  | 1108570 | 8.5x70 mm  |
| 1104025 | 4.0x25 mm | 1106025 | 6.0x25 mm | 1107580 | 7.5x80 mm  | 1108580 | 8.5x80 mm  |
| 1104030 | 4.0x30 mm | 1106030 | 6.0x30 mm | 1107590 | 7.5x90 mm  | 1108590 | 8.5x90 mm  |
| 1104035 | 4.0x35 mm | 1106035 | 6.0x35 mm | 1108025 | 8.0x25 mm  | 1108510 | 8.5x100 mm |
| 1104040 | 4.0x40 mm | 1106040 | 6.0x40 mm | 1108030 | 8.0x30 mm  |         |            |
| 1104045 | 4.0x45 mm | 1106045 | 6.0x45 mm | 1108035 | 8.0x35 mm  |         |            |
| 1104050 | 4.0x50 mm | 1106050 | 6.0x50 mm | 1108040 | 8.0x40 mm  |         |            |
| 1104055 | 4.0x55 mm | 1106055 | 6.0x55 mm | 1108045 | 8.0x45 mm  |         |            |
| 1104520 | 4.5x20 mm | 1106060 | 6.0x60 mm | 1108050 | 8.0x50 mm  |         |            |
| 1104525 | 4.5x25 mm | 1106070 | 6.0x70 mm | 1108055 | 8.0x55 mm  |         |            |
| 1104530 | 4.5x30 mm | 1106525 | 6.5x25 mm | 1108060 | 8.0x60 mm  |         |            |
| 1104535 | 4.5x35 mm | 1106530 | 6.5x30 mm | 1108070 | 8.0x70 mm  |         |            |
| 1104540 | 4.5x40 mm | 1106535 | 6.5x35 mm | 1108080 | 8.0x80 mm  |         |            |
| 1104545 | 4.5x45 mm | 1106540 | 6.5x40 mm | 1108090 | 8.0x90 mm  |         |            |
| 1104550 | 4.5x50 mm | 1106545 | 6.5x45 mm | 1108010 | 8.0x100 mm |         |            |
| 1104555 | 4.5x55 mm | 1106550 | 6.5x50 mm |         |            |         |            |
| 1105020 | 5.0x20 mm | 1106555 | 6.5x55 mm |         |            |         |            |
| 1105025 | 5.0x25 mm | 1106560 | 6.5x60 mm |         |            |         |            |
| 1105030 | 5.0x30 mm | 1106570 | 6.5x70 mm |         |            |         |            |
| 1105035 | 5.0x35 mm | 1107025 | 7.0x25 mm |         |            |         |            |
| 1105040 | 5.0x40 mm | 1107030 | 7.0x30 mm |         |            |         |            |
| 1105045 | 5.0x45 mm | 1107035 | 7.0x35 mm |         |            |         |            |
| 1105050 | 5.0x50 mm | 1107040 | 7.0x40 mm |         |            |         |            |
| 1105055 | 5.0x55 mm | 1107045 | 7.0x45 mm |         |            |         |            |
| 1105060 | 5.0x60 mm | 1107050 | 7.0x50 mm |         |            |         |            |
|         |           | 1107055 | 7.0x55 mm |         |            |         |            |
|         |           | 1107060 | 7.0x60 mm |         |            |         |            |
|         |           | 1107070 | 7.0x70 mm |         |            |         |            |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:
  1. Cannulated screw allowing cement injection
  2. Double-threaded Screw providing a tighter grip
  3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



## Posterior Spinal Screw System

### 4U SPINAL POSTERIOR MULTIFUNCTIONAL PEDICLE SCREW



| Code    | Size      | Code    | Size       | Code    | Size       |
|---------|-----------|---------|------------|---------|------------|
| 1203525 | 3.5x25 mm | 1206530 | 6.5x30 mm  | 1208530 | 8.5x30 mm  |
| 1203530 | 3.5x30 mm | 1206535 | 6.5x35 mm  | 1208535 | 8.5x35 mm  |
| 1203535 | 3.5x35 mm | 1206540 | 6.5x40 mm  | 1208540 | 8.5x40 mm  |
| 1203540 | 3.5x40 mm | 1206545 | 6.5x45 mm  | 1208545 | 8.5x45 mm  |
| 1203545 | 3.5x45 mm | 1206550 | 6.5x50 mm  | 1208550 | 8.5x50 mm  |
| 1204025 | 4.0x25 mm | 1206555 | 6.5x55 mm  | 1208555 | 8.5x55 mm  |
| 1204030 | 4.0x30 mm | 1206560 | 6.5x60 mm  | 1208560 | 8.5x60 mm  |
| 1204035 | 4.0x35 mm | 1207030 | 7.0x30 mm  | 1208565 | 8.5x65 mm  |
| 1204040 | 4.0x40 mm | 1207035 | 7.0x35 mm  | 1208570 | 8.5x70 mm  |
| 1204045 | 4.0x45 mm | 1207040 | 7.0x40 mm  | 1208580 | 8.5x80 mm  |
| 1204525 | 4.5x25 mm | 1207045 | 7.0x45 mm  | 1208590 | 8.5x90 mm  |
| 1204530 | 4.5x30 mm | 1207050 | 7.0x50 mm  | 1208510 | 8.5x100 mm |
| 1204535 | 4.5x35 mm | 1207055 | 7.0x55 mm  |         |            |
| 1204540 | 4.5x40 mm | 1207060 | 7.0x60 mm  |         |            |
| 1204545 | 4.5x45 mm | 1207530 | 7.5x30 mm  |         |            |
| 1205030 | 5.0x30 mm | 1207535 | 7.5x35 mm  |         |            |
| 1205035 | 5.0x35 mm | 1207540 | 7.5x40 mm  |         |            |
| 1205040 | 5.0x40 mm | 1207545 | 7.5x45 mm  |         |            |
| 1205045 | 5.0x45 mm | 1207550 | 7.5x50 mm  |         |            |
| 1205050 | 5.0x50 mm | 1207555 | 7.5x55 mm  |         |            |
| 1205055 | 5.0x55 mm | 1207560 | 7.5x60 mm  |         |            |
| 1205525 | 5.5x25 mm | 1208030 | 8.0x30 mm  |         |            |
| 1205530 | 5.5x30 mm | 1208035 | 8.0x35 mm  |         |            |
| 1205535 | 5.5x35 mm | 1208040 | 8.0x40 mm  |         |            |
| 1205540 | 5.5x40 mm | 1208045 | 8.0x45 mm  |         |            |
| 1205545 | 5.5x45 mm | 1208050 | 8.0x50 mm  |         |            |
| 1205550 | 5.5x50 mm | 1208055 | 8.0x55 mm  |         |            |
| 1205555 | 5.5x55 mm | 1208060 | 8.0x60 mm  |         |            |
| 1205560 | 5.5x60 mm | 1208065 | 8.0x65 mm  |         |            |
| 1206030 | 6.0x30 mm | 1208070 | 8.0x70 mm  |         |            |
| 1206035 | 6.0x35 mm | 1208080 | 8.0x80 mm  |         |            |
| 1206040 | 6.0x40 mm | 1208090 | 8.0x90 mm  |         |            |
| 1206045 | 6.0x45 mm | 1208010 | 8.0x100 mm |         |            |
| 1206050 | 6.0x50 mm |         |            |         |            |
| 1206055 | 6.0x55 mm |         |            |         |            |
| 1206060 | 6.0x60 mm |         |            |         |            |

# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:

1. Cannulated screw allowing cement injection
2. Double-threaded Screw providing a tighter grip
3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine

## Posterior Spinal Screw System

### 4U SPINAL PEDICLE POSTERIOR MONOAXIAL SCREW


| Code    | Size      | Code    | Size      |
|---------|-----------|---------|-----------|
| 1303525 | 3.5x25 mm | 1306035 | 6.0x35 mm |
| 1303530 | 3.5x30 mm | 1306040 | 6.0x40 mm |
| 1303535 | 3.5x35 mm | 1306045 | 6.0x45 mm |
| 1303540 | 3.5x40 mm | 1306050 | 6.0x50 mm |
| 1303545 | 3.5x45 mm | 1306055 | 6.0x55 mm |
| 1304025 | 4.0x25 mm | 1306530 | 6.5x30 mm |
| 1304030 | 4.0x30 mm | 1306535 | 6.5x35 mm |
| 1304035 | 4.0x35 mm | 1306540 | 6.5x40 mm |
| 1304040 | 4.0x40 mm | 1306545 | 6.5x45 mm |
| 1304045 | 4.0x45 mm | 1306550 | 6.5x50 mm |
| 1304525 | 4.5x25 mm | 1306555 | 6.5x55 mm |
| 1304530 | 4.5x30 mm | 1307035 | 7.0x35 mm |
| 1304535 | 4.5x35 mm | 1307040 | 7.0x40 mm |
| 1304540 | 4.5x40 mm | 1307045 | 7.0x45 mm |
| 1304545 | 4.5x45 mm | 1307050 | 7.0x50 mm |
| 1305030 | 5.0x30 mm | 1307055 | 7.0x55 mm |
| 1305035 | 5.0x35 mm | 1307535 | 7.5x35 mm |
| 1305040 | 5.0x40 mm | 1307540 | 7.5x40 mm |
| 1305045 | 5.0x45 mm | 1307545 | 7.5x45 mm |
| 1305050 | 5.0x50 mm | 1307550 | 7.5x50 mm |
| 1305530 | 5.5x30 mm | 1307555 | 7.5x55 mm |
| 1305535 | 5.5x35 mm | 1308035 | 8.0x35 mm |
| 1305540 | 5.5x40 mm | 1308040 | 8.0x40 mm |
| 1305545 | 5.5x45 mm | 1308045 | 8.0x45 mm |
| 1305550 | 5.5x50 mm | 1308050 | 8.0x50 mm |
| 1305555 | 5.5x55 mm | 1308055 | 8.0x55 mm |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:
  1. Cannulated screw allowing cement injection
  2. Double-threaded Screw providing a tighter grip
  3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



## Posterior Spinal Screw System

### 4U SPINAL PEDICLE POSTERIOR MONOAXIAL SCREW

#### CANNULATED

| Code    | Size      |
|---------|-----------|
| 1315535 | 5.5x35 mm |
| 1315540 | 5.5x40 mm |
| 1315545 | 5.5x45 mm |
| 1315550 | 5.5x50 mm |
| 1315555 | 5.5x55 mm |
| 1316530 | 6.5x30 mm |
| 1316535 | 6.5x35 mm |
| 1316540 | 6.5x40 mm |
| 1316545 | 6.5x45 mm |
| 1316550 | 6.5x50 mm |
| 1316555 | 6.5x55 mm |
| 1317530 | 7.5x35 mm |
| 1317535 | 7.5x35 mm |
| 1317540 | 7.5x40 mm |
| 1317545 | 7.5x45 mm |
| 1317550 | 7.5x50 mm |
| 1317555 | 7.5x55 mm |



#### CANNULATED-CEMENT


| Code    | Size      |
|---------|-----------|
| 1325535 | 5.5x35 mm |
| 1325540 | 5.5x40 mm |
| 1325545 | 5.5x45 mm |
| 1325550 | 5.5x50 mm |
| 1325555 | 5.5x55 mm |
| 1326530 | 6.5x30 mm |
| 1326535 | 6.5x35 mm |
| 1326540 | 6.5x40 mm |
| 1326545 | 6.5x45 mm |
| 1326550 | 6.5x50 mm |
| 1326555 | 6.5x55 mm |
| 1327530 | 7.5x35 mm |
| 1327535 | 7.5x35 mm |
| 1327540 | 7.5x40 mm |
| 1327545 | 7.5x45 mm |
| 1327550 | 7.5x50 mm |
| 1327555 | 7.5x55 mm |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:
  1. Cannulated screw allowing cement injection
  2. Double-threaded Screw providing a tighter grip
  3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



## Posterior Spinal Screw System

### 4U SPINAL PEDICLE POSTERIOR POLYAXIAL CANNULATED SCREW


| Code    | Size      | Code    | Size      | Code    | Size      |
|---------|-----------|---------|-----------|---------|-----------|
| 1604030 | 4.0x30 mm | 1605530 | 5.5x30 mm | 1607030 | 7.0x30 mm |
| 1604035 | 4.0x35 mm | 1605535 | 5.5x35 mm | 1607035 | 7.0x35 mm |
| 1604040 | 4.0x40 mm | 1605540 | 5.5x40 mm | 1607040 | 7.0x40 mm |
| 1604045 | 4.0x45 mm | 1605545 | 5.5x45 mm | 1607045 | 7.0x45 mm |
| 1604050 | 4.0x50 mm | 1605550 | 5.5x50 mm | 1607050 | 7.0x50 mm |
| 1604055 | 4.0x55 mm | 1605555 | 5.5x55 mm | 1607055 | 7.0x55 mm |
| 1604060 | 4.0x60 mm | 1605560 | 5.5x60 mm | 1607060 | 7.0x60 mm |
| 1604065 | 4.0x65 mm | 1605565 | 5.5x65 mm | 1607065 | 7.0x65 mm |
| 1604070 | 4.0x70 mm | 1605570 | 5.5x70 mm | 1607070 | 7.0x70 mm |
| 1604530 | 4.5x30 mm | 1606030 | 6.0x30 mm | 1607080 | 7.0x80 mm |
| 1604535 | 4.5x35 mm | 1606035 | 6.0x35 mm | 1607530 | 7.5x30 mm |
| 1604540 | 4.5x40 mm | 1606040 | 6.0x40 mm | 1607535 | 7.5x35 mm |
| 1604545 | 4.5x45 mm | 1606045 | 6.0x45 mm | 1607540 | 7.5x40 mm |
| 1604550 | 4.5x50 mm | 1606050 | 6.0x50 mm | 1607545 | 7.5x45 mm |
| 1604555 | 4.5x55 mm | 1606055 | 6.0x55 mm | 1607550 | 7.5x50 mm |
| 1604560 | 4.5x60 mm | 1606060 | 6.0x60 mm | 1607555 | 7.5x55 mm |
| 1604565 | 4.5x65 mm | 1606065 | 6.0x65 mm | 1607560 | 7.5x60 mm |
| 1604570 | 4.5x70 mm | 1606070 | 6.0x70 mm | 1607565 | 7.5x65 mm |
| 1604580 | 4.5x80 mm | 1606080 | 6.0x80 mm | 1607570 | 7.5x70 mm |
| 1605030 | 5.0x30 mm | 1606530 | 6.5x30 mm | 1607580 | 7.5x80 mm |
| 1605035 | 5.0x35 mm | 1606535 | 6.5x35 mm | 1608035 | 8.0x35 mm |
| 1605040 | 5.0x40 mm | 1606540 | 6.5x40 mm | 1608040 | 8.0x40 mm |
| 1605045 | 5.0x45 mm | 1606545 | 6.5x45 mm | 1608045 | 8.0x45 mm |
| 1605050 | 5.0x50 mm | 1606550 | 6.5x50 mm | 1608050 | 8.0x50 mm |
| 1605055 | 5.0x55 mm | 1606555 | 6.5x55 mm | 1608055 | 8.0x55 mm |
| 1605060 | 5.0x60 mm | 1606560 | 6.5x60 mm | 1608060 | 8.0x60 mm |
| 1605065 | 5.0x65 mm | 1606565 | 6.5x65 mm | 1608065 | 8.0x65 mm |
| 1605070 | 5.0x70 mm | 1606570 | 6.5x70 mm | 1608070 | 8.0x70 mm |
| 1605080 | 5.0x80 mm | 1606580 | 6.5x80 mm | 1608080 | 8.0x80 mm |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:
  1. Cannulated screw allowing cement injection
  2. Double-threaded Screw providing a tighter grip
  3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



## Posterior Spinal Screw System

### 4U SPINAL PEDICLE POSTERIOR POLYAXIAL CANNULATED-CEMENT SCREW

| Code    | Size      | Code    | Size      | Code    | Size      |
|---------|-----------|---------|-----------|---------|-----------|
| 1704030 | 4.0x30 mm | 1705530 | 5.5x30 mm | 1707030 | 7.0x30 mm |
| 1704035 | 4.0x35 mm | 1705535 | 5.5x35 mm | 1707035 | 7.0x35 mm |
| 1704040 | 4.0x40 mm | 1705540 | 5.5x40 mm | 1707040 | 7.0x40 mm |
| 1704045 | 4.0x45 mm | 1705545 | 5.5x45 mm | 1707045 | 7.0x45 mm |
| 1704050 | 4.0x50 mm | 1705550 | 5.5x50 mm | 1707050 | 7.0x50 mm |
| 1704055 | 4.0x55 mm | 1705555 | 5.5x55 mm | 1707055 | 7.0x55 mm |
| 1704060 | 4.0x60 mm | 1705560 | 5.5x60 mm | 1707060 | 7.0x60 mm |
| 1704065 | 4.0x65 mm | 1705565 | 5.5x65 mm | 1707065 | 7.0x65 mm |
| 1704070 | 4.0x70 mm | 1705570 | 5.5x70 mm | 1707070 | 7.0x70 mm |
| 1704530 | 4.5x30 mm | 1706030 | 6.0x30 mm | 1707080 | 7.0x80 mm |
| 1704535 | 4.5x35 mm | 1706035 | 6.0x35 mm | 1707530 | 7.5x30 mm |
| 1704540 | 4.5x40 mm | 1706040 | 6.0x40 mm | 1707535 | 7.5x35 mm |
| 1704545 | 4.5x45 mm | 1706045 | 6.0x45 mm | 1707540 | 7.5x40 mm |
| 1704550 | 4.5x50 mm | 1706050 | 6.0x50 mm | 1707545 | 7.5x45 mm |
| 1704555 | 4.5x55 mm | 1706055 | 6.0x55 mm | 1707550 | 7.5x50 mm |
| 1704560 | 4.5x60 mm | 1706060 | 6.0x60 mm | 1707555 | 7.5x55 mm |
| 1704565 | 4.5x65 mm | 1706065 | 6.0x65 mm | 1707560 | 7.5x60 mm |
| 1704570 | 4.5x70 mm | 1706070 | 6.0x70 mm | 1707565 | 7.5x65 mm |
| 1704580 | 4.5x80 mm | 1706080 | 6.0x80 mm | 1707570 | 7.5x70 mm |
| 1705030 | 5.0x30 mm | 1706530 | 6.5x30 mm | 1707580 | 7.5x80 mm |
| 1705035 | 5.0x35 mm | 1706535 | 6.5x35 mm | 1708035 | 8.0x35 mm |
| 1705040 | 5.0x40 mm | 1706540 | 6.5x40 mm | 1708040 | 8.0x40 mm |
| 1705045 | 5.0x45 mm | 1706545 | 6.5x45 mm | 1708045 | 8.0x45 mm |
| 1705050 | 5.0x50 mm | 1706550 | 6.5x50 mm | 1708050 | 8.0x50 mm |
| 1705055 | 5.0x55 mm | 1706555 | 6.5x55 mm | 1708055 | 8.0x55 mm |
| 1705060 | 5.0x60 mm | 1706560 | 6.5x60 mm | 1708060 | 8.0x60 mm |
| 1705065 | 5.0x65 mm | 1706565 | 6.5x65 mm | 1708065 | 8.0x65 mm |
| 1705070 | 5.0x70 mm | 1706570 | 6.5x70 mm | 1708070 | 8.0x70 mm |
| 1705080 | 5.0x80 mm | 1706580 | 6.5x80 mm | 1708080 | 8.0x80 mm |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:

1. Cannulated screw allowing cement injection
2. Double-threaded Screw providing a tighter grip
3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



**4U Spinal Posterior Polyaxial Spondylolisthesis Screw**

## 4U SPINAL POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW

| Code    | Size      | Code    | Size      |
|---------|-----------|---------|-----------|
| 1404030 | 4.0x30 mm | 1406535 | 6.5x35 mm |
| 1404035 | 4.0x35 mm | 1406540 | 6.5x40 mm |
| 1404040 | 4.0x40 mm | 1406545 | 6.5x45 mm |
| 1404045 | 4.0x45 mm | 1406550 | 6.5x50 mm |
| 1404050 | 4.0x50 mm | 1406555 | 6.5x55 mm |
| 1404530 | 4.5x30 mm | 1406560 | 6.5x60 mm |
| 1404535 | 4.5x35 mm | 1406565 | 6.5x65 mm |
| 1404540 | 4.5x40 mm | 1406570 | 6.5x70 mm |
| 1404545 | 4.5x45 mm | 1407035 | 7.0x35 mm |
| 1404550 | 4.5x50 mm | 1407040 | 7.0x40 mm |
| 1405035 | 5.0x35 mm | 1407045 | 7.0x45 mm |
| 1405040 | 5.0x40 mm | 1407050 | 7.0x50 mm |
| 1405045 | 5.0x45 mm | 1407055 | 7.0x55 mm |
| 1405050 | 5.0x50 mm | 1407060 | 7.0x60 mm |
| 1405055 | 5.0x55 mm | 1407065 | 7.0x65 mm |
| 1405060 | 5.0x60 mm | 1407070 | 7.0x70 mm |
| 1405525 | 5.5x25 mm | 1407535 | 7.5x35 mm |
| 1405530 | 5.5x30 mm | 1407540 | 7.5x40 mm |
| 1405535 | 5.5x35 mm | 1407545 | 7.5x45 mm |
| 1405540 | 5.5x40 mm | 1407550 | 7.5x50 mm |
| 1405545 | 5.5x45 mm | 1407555 | 7.5x55 mm |
| 1405550 | 5.5x50 mm | 1407560 | 7.5x60 mm |
| 1405555 | 5.5x55 mm | 1407565 | 7.5x65 mm |
| 1405560 | 5.5x60 mm | 1407570 | 7.5x70 mm |
| 1405565 | 5.5x65 mm | 1407580 | 7.5x80 mm |
| 1406035 | 6.0x35 mm | 1408035 | 7.0x35 mm |
| 1406040 | 6.0x40 mm | 1408040 | 7.0x40 mm |
| 1406045 | 6.0x45 mm | 1408045 | 7.0x45 mm |
| 1406050 | 6.0x50 mm | 1408050 | 7.0x50 mm |
| 1406055 | 6.0x55 mm | 1408055 | 7.0x55 mm |
| 1406060 | 6.0x60 mm | 1408060 | 7.0x60 mm |
| 1406065 | 6.0x65 mm | 1408065 | 7.0x65 mm |
|         |           | 1408070 | 7.0x70 mm |
|         |           | 1408080 | 7.0x80 mm |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:

1. Cannulated screw allowing cement injection
2. Double-threaded Screw providing a tighter grip
3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



**4U Spinal Posterior Polyaxial Spondylolisthesis Screw**

## 4U SPINAL POSTERIOR POLYAXIAL SPONDYLOLISTHESIS SCREW

### MULTIFUNCTIONAL

### CANNULATED

### CANNULATED-CEMENT

#### Code Size

|         |           |
|---------|-----------|
| 1434530 | 4.5x30 mm |
| 1434535 | 4.5x35 mm |
| 1434540 | 4.5x40 mm |
| 1434545 | 4.5x45 mm |
| 1435535 | 5.5x35 mm |
| 1435540 | 5.5x40 mm |
| 1435545 | 5.5x45 mm |
| 1435550 | 5.5x50 mm |
| 1436035 | 6.0x35 mm |
| 1436040 | 6.0x40 mm |
| 1436045 | 6.0x45 mm |
| 1436050 | 6.0x50 mm |
| 1436535 | 6.5x35 mm |
| 1436540 | 6.5x40 mm |
| 1436545 | 6.5x45 mm |
| 1436550 | 6.5x50 mm |
| 1437035 | 7.0x35 mm |
| 1437040 | 7.0x40 mm |
| 1437045 | 7.0x45 mm |
| 1437050 | 7.0x50 mm |
| 1437055 | 7.0x55 mm |
| 1437540 | 7.5x40 mm |
| 1437545 | 7.5x45 mm |
| 1437550 | 7.5x50 mm |
| 1437555 | 7.5x55 mm |



#### Code Size

|         |           |
|---------|-----------|
| 1414530 | 4.5x30 mm |
| 1414535 | 4.5x35 mm |
| 1414540 | 4.5x40 mm |
| 1414545 | 4.5x45 mm |
| 1415535 | 5.5x35 mm |
| 1415540 | 5.5x40 mm |
| 1415545 | 5.5x45 mm |
| 1415550 | 5.5x50 mm |
| 1416035 | 6.0x35 mm |
| 1416040 | 6.0x40 mm |
| 1416045 | 6.0x45 mm |
| 1416050 | 6.0x50 mm |
| 1416535 | 6.5x35 mm |
| 1416540 | 6.5x40 mm |
| 1416545 | 6.5x45 mm |
| 1416550 | 6.5x50 mm |
| 1417035 | 7.0x35 mm |
| 1417040 | 7.0x40 mm |
| 1417045 | 7.0x45 mm |
| 1417050 | 7.0x50 mm |
| 1417055 | 7.0x55 mm |
| 1417540 | 7.5x40 mm |
| 1417545 | 7.5x45 mm |
| 1417550 | 7.5x50 mm |
| 1417555 | 7.5x55 mm |



#### Code Size

|         |           |
|---------|-----------|
| 1424530 | 4.5x30 mm |
| 1424535 | 4.5x35 mm |
| 1424540 | 4.5x40 mm |
| 1424545 | 4.5x45 mm |
| 1425535 | 5.5x35 mm |
| 1425540 | 5.5x40 mm |
| 1425545 | 5.5x45 mm |
| 1425550 | 5.5x50 mm |
| 1426035 | 6.0x35 mm |
| 1426040 | 6.0x40 mm |
| 1426045 | 6.0x45 mm |
| 1426050 | 6.0x50 mm |
| 1426535 | 6.5x35 mm |
| 1426540 | 6.5x40 mm |
| 1426545 | 6.5x45 mm |
| 1426550 | 6.5x50 mm |
| 1427035 | 7.0x35 mm |
| 1427040 | 7.0x40 mm |
| 1427045 | 7.0x45 mm |
| 1427050 | 7.0x50 mm |
| 1427055 | 7.0x55 mm |
| 1427540 | 7.5x40 mm |
| 1427545 | 7.5x45 mm |
| 1427550 | 7.5x50 mm |
| 1427555 | 7.5x55 mm |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial Cannulated screws, various types of designs for different applications:

1. Cannulated screw allowing cement injection
2. Double-threaded Screw providing a tighter grip
3. Reduction Screw which is specifically designed to address spondylolisthesis and other deformities of the spine



**4U Spinal Posterior Monoaxial Spondylolisthesis Screw**

## 4U SPINAL POSTERIOR MONOAXIAL SPONDYLOLISTHESIS SCREW

### STANDART

### CANNULATED

### CANNULATED-CEMENT

#### Code

#### Size

#### Code

#### Size

#### Code

#### Size



1504530 4.5x30 mm  
1504535 4.5x35 mm  
1504540 4.5x40 mm  
1504545 4.5x45 mm  
1505535 5.5x35 mm  
1505540 5.5x40 mm  
1505545 5.5x45 mm  
1505550 5.5x50 mm  
1506035 6.0x35 mm  
1506040 6.0x40 mm  
1506045 6.0x45 mm  
1506050 6.0x50 mm  
1506535 6.5x35 mm  
1506540 6.5x40 mm  
1506545 6.5x45 mm  
1506550 6.5x50 mm  
1507035 7.0x35 mm  
1507040 7.0x40 mm  
1507045 7.0x45 mm  
1507050 7.0x50 mm  
1507055 7.0x55 mm  
1507540 7.5x40 mm  
1507545 7.5x45 mm  
1507550 7.5x50 mm  
1507555 7.5x55 mm



1514530 4.5x30 mm  
1514535 4.5x35 mm  
1514540 4.5x40 mm  
1514545 4.5x45 mm  
1515535 5.5x35 mm  
1515540 5.5x40 mm  
1515545 5.5x45 mm  
1515550 5.5x50 mm  
1516035 6.0x35 mm  
1516040 6.0x40 mm  
1516045 6.0x45 mm  
1516050 6.0x50 mm  
1516535 6.5x35 mm  
1516540 6.5x40 mm  
1516545 6.5x45 mm  
1516550 6.5x50 mm  
1517035 7.0x35 mm  
1517040 7.0x40 mm  
1517045 7.0x45 mm  
1517050 7.0x50 mm  
1517055 7.0x55 mm  
1517540 7.5x40 mm  
1517545 7.5x45 mm  
1517550 7.5x50 mm  
1517555 7.5x55 mm



1524530 4.5x30 mm  
1524535 4.5x35 mm  
1524540 4.5x40 mm  
1524545 4.5x45 mm  
1525535 5.5x35 mm  
1525540 5.5x40 mm  
1525545 5.5x45 mm  
1525550 5.5x50 mm  
1526035 6.0x35 mm  
1526040 6.0x40 mm  
1526045 6.0x45 mm  
1526050 6.0x50 mm  
1526535 6.5x35 mm  
1526540 6.5x40 mm  
1526545 6.5x45 mm  
1526550 6.5x50 mm  
1527035 7.0x35 mm  
1527040 7.0x40 mm  
1527045 7.0x45 mm  
1527050 7.0x50 mm  
1527055 7.0x55 mm  
1527540 7.5x40 mm  
1527545 7.5x45 mm  
1527550 7.5x50 mm  
1527555 7.5x55 mm

# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial, Cannulated, Cannulated Cemented screws, various types of designs for different applications
- Ability to accept  $\Phi 5.5$ mm rods
- Flexible move in 360° direction
- Curled spherical head provides better angular stability
- Multi-type extended tap design
- Double-lead thread pattern for faster insert on and increased pullout strength
- Serrated saddle locking block prevents the loosening of screw.



## 4U Spinal Pedicle Posterior Polyaxial Screw Long Arm

### 4U SPINAL PEDICLE POSTERIOR POLYAXIAL SCREW LONG ARM

#### STANDART

#### MULTIFUNCTIONAL

| Code    | Size      | Code    | Size      | Code    | Size      | Code    | Size      |
|---------|-----------|---------|-----------|---------|-----------|---------|-----------|
| 1675530 | 5.5X30 MM | 1677030 | 7.0X30 MM | 1685530 | 5.5X30 MM | 1687030 | 7.0X30 MM |
| 1675535 | 5.5X35 MM | 1677035 | 7.0X35 MM | 1685535 | 5.5X35 MM | 1687035 | 7.0X35 MM |
| 1675540 | 5.5X40 MM | 1677040 | 7.0X40 MM | 1685540 | 5.5X40 MM | 1687040 | 7.0X40 MM |
| 1675545 | 5.5X45 MM | 1677045 | 7.0X45 MM | 1685545 | 5.5X45 MM | 1687045 | 7.0X45 MM |
| 1675550 | 5.5X50 MM | 1677050 | 7.0X50 MM | 1685550 | 5.5X50 MM | 1687050 | 7.0X50 MM |
| 1675555 | 5.5X55 MM | 1677055 | 7.0X55 MM | 1685555 | 5.5X55 MM | 1687055 | 7.0X55 MM |
| 1676030 | 6.0X30 MM | 1677530 | 7.5X30 MM | 1686030 | 6.0X30 MM | 1687530 | 7.5X30 MM |
| 1676035 | 6.0X35 MM | 1677535 | 7.5X35 MM | 1686035 | 6.0X35 MM | 1687535 | 7.5X35 MM |
| 1676040 | 6.0X40 MM | 1677540 | 7.5X40 MM | 1686040 | 6.0X40 MM | 1687540 | 7.5X40 MM |
| 1676045 | 6.0X45 MM | 1677545 | 7.5X45 MM | 1686045 | 6.0X45 MM | 1687545 | 7.5X45 MM |
| 1676050 | 6.0X50 MM | 1677550 | 7.5X50 MM | 1686050 | 6.0X50 MM | 1687550 | 7.5X50 MM |
| 1676055 | 6.0X55 MM | 1677555 | 7.5X55 MM | 1686055 | 6.0X55 MM | 1687555 | 7.5X55 MM |
| 1676530 | 6.5X30 MM | 1678030 | 8.0X30 MM | 1686530 | 6.5X30 MM | 1688030 | 8.0X30 MM |
| 1676535 | 6.5X35 MM | 1678035 | 8.0X35 MM | 1686535 | 6.5X35 MM | 1688035 | 8.0X35 MM |
| 1676540 | 6.5X40 MM | 1678040 | 8.0X40 MM | 1686540 | 6.5X40 MM | 1688040 | 8.0X40 MM |
| 1676545 | 6.5X45 MM | 1678045 | 8.0X45 MM | 1686545 | 6.5X45 MM | 1688045 | 8.0X45 MM |
| 1676550 | 6.5X50 MM | 1678050 | 8.0X50 MM | 1686550 | 6.5X50 MM | 1688050 | 8.0X50 MM |
| 1676555 | 6.5X55 MM | 1678055 | 8.0X55 MM | 1686555 | 6.5X55 MM | 1688055 | 8.0X55 MM |



# Posterior Spinal Screw System

## Features & Benefits

- Ti6Al4V-ELI raw material
- Low profile screws
- 57° Angulation in polyaxial screws
- Cylindrical and self-tapping feature
- Easy screw fixation
- Flexibility for screw and rod replacement
- Monoaxial, Polyaxial, Cannulated, Cannulated Cemented screws, various types of designs for different applications
- Ability to accept Ø5.5mm rods
- Flexible move in 360° direction
- Curled spherical head provides better angular stability
- Multy-type extended tap design
- Double-lead thread pattern for faster insert on and increased pullout strength
- Serrated saddle locking block prevents the loosening of screw.

## 4U Spinal Pedicle Posterior Polyaxial Screw Long Arm

### 4U SPINAL PEDICLE POSTERIOR POLYAXIAL SCREW LONG ARM

#### CANNULATED

#### CANNULATED-CEMENT


| Code    | Size      | Code    | Size      | Code    | Size      | Code    | Size      |
|---------|-----------|---------|-----------|---------|-----------|---------|-----------|
| 1655530 | 5.5X30 MM | 1657030 | 7.0X30 MM | 1665530 | 5.5X30 MM | 1667030 | 7.0X30 MM |
| 1655535 | 5.5X35 MM | 1657035 | 7.0X35 MM | 1665535 | 5.5X35 MM | 1667035 | 7.0X35 MM |
| 1655540 | 5.5X40 MM | 1657040 | 7.0X40 MM | 1665540 | 5.5X40 MM | 1667040 | 7.0X40 MM |
| 1655545 | 5.5X45 MM | 1657045 | 7.0X45 MM | 1665545 | 5.5X45 MM | 1667045 | 7.0X45 MM |
| 1655550 | 5.5X50 MM | 1657050 | 7.0X50 MM | 1665550 | 5.5X50 MM | 1667050 | 7.0X50 MM |
| 1655555 | 5.5X55 MM | 1657055 | 7.0X55 MM | 1665555 | 5.5X55 MM | 1667055 | 7.0X55 MM |
| 1656030 | 6.0X30 MM | 1657530 | 7.5X30 MM | 1666030 | 6.0X30 MM | 1667530 | 7.5X30 MM |
| 1656035 | 6.0X35 MM | 1657535 | 7.5X35 MM | 1666035 | 6.0X35 MM | 1667535 | 7.5X35 MM |
| 1656040 | 6.0X40 MM | 1657540 | 7.5X40 MM | 1666040 | 6.0X40 MM | 1667540 | 7.5X40 MM |
| 1656045 | 6.0X45 MM | 1657545 | 7.5X45 MM | 1666045 | 6.0X45 MM | 1667545 | 7.5X45 MM |
| 1656050 | 6.0X50 MM | 1657550 | 7.5X50 MM | 1666050 | 6.0X50 MM | 1667550 | 7.5X50 MM |
| 1656055 | 6.0X55 MM | 1657555 | 7.5X55 MM | 1666055 | 6.0X55 MM | 1667555 | 7.5X55 MM |
| 1656530 | 6.5X30 MM | 1658030 | 8.0X30 MM | 1666530 | 6.5X30 MM | 1668030 | 8.0X30 MM |
| 1656535 | 6.5X35 MM | 1658035 | 8.0X35 MM | 1666535 | 6.5X35 MM | 1668035 | 8.0X35 MM |
| 1656540 | 6.5X40 MM | 1658040 | 8.0X40 MM | 1666540 | 6.5X40 MM | 1668040 | 8.0X40 MM |
| 1656545 | 6.5X45 MM | 1658045 | 8.0X45 MM | 1666545 | 6.5X45 MM | 1668045 | 8.0X45 MM |
| 1656550 | 6.5X50 MM | 1658050 | 8.0X50 MM | 1666550 | 6.5X50 MM | 1668050 | 8.0X50 MM |
| 1656555 | 6.5X55 MM | 1658055 | 8.0X55 MM | 1666555 | 6.5X55 MM | 1668055 | 8.0X55 MM |



# Posterior Spinal Rods

## Features & Benefits

- Ti6Al4V-ELI raw material
- Easy screw fixation
- Flexibility for rod replacement
- Appropriate straight & curved rod in a variety of lengths
- Alternatively rods may be cut to length using the cutter
- The straight rod can be curved the desired lordosis using the Rod Bender
- 5.5 and 6.0 option



**Posterior Spinal Rods, 5.5**

### Spinal Rod 5.5

| Code    | Size       | Code    | Size       |
|---------|------------|---------|------------|
| 2300540 | 5.5x40 mm  | 2305270 | 5.5x270 mm |
| 2300550 | 5.5x50 mm  | 2305280 | 5.5x280 mm |
| 2300560 | 5.5x60 mm  | 2305290 | 5.5x290 mm |
| 2300570 | 5.5x70 mm  | 2305300 | 5.5x300 mm |
| 2300580 | 5.5x80 mm  | 2305310 | 5.5x310 mm |
| 2300590 | 5.5x90 mm  | 2305320 | 5.5x320 mm |
| 2305100 | 5.5x100 mm | 2305340 | 5.5x340 mm |
| 2305110 | 5.5x110 mm | 2305350 | 5.5x350 mm |
| 2305120 | 5.5x120 mm | 2305360 | 5.5x360 mm |
| 2305130 | 5.5x130 mm | 2305380 | 5.5x380 mm |
| 2305140 | 5.5x140 mm | 2305400 | 5.5x400 mm |
| 2305150 | 5.5x150 mm | 2305420 | 5.5x420 mm |
| 2305160 | 5.5x160 mm | 2305440 | 5.5x440 mm |
| 2305170 | 5.5x170 mm | 2305450 | 5.5x450 mm |
| 2305180 | 5.5x180 mm | 2305460 | 5.5x460 mm |
| 2305190 | 5.5x190 mm | 2305480 | 5.5x480 mm |
| 2305200 | 5.5x200 mm | 2305500 | 5.5x500 mm |
| 2305210 | 5.5x210 mm | 2305520 | 5.5x520 mm |
| 2305220 | 5.5x220 mm | 2305540 | 5.5x540 mm |
| 2305230 | 5.5x230 mm | 2305550 | 5.5x550 mm |
| 2305240 | 5.5x240 mm | 2305560 | 5.5x560 mm |
| 2305250 | 5.5x250 mm | 2305580 | 5.5x580 mm |
| 2305260 | 5.5x260 mm | 2305600 | 5.5x600 mm |


### 4U SPINAL HIGH FLEX ROD 5.5

| Code    | Size       |
|---------|------------|
| 2155035 | 5.5x35 mm  |
| 2155040 | 5.5x40 mm  |
| 2155045 | 5.5x45 mm  |
| 2155050 | 5.5x50 mm  |
| 2155055 | 5.5x55 mm  |
| 2155060 | 5.5x60 mm  |
| 2155065 | 5.5x65 mm  |
| 2155070 | 5.5x70 mm  |
| 2155080 | 5.5x80 mm  |
| 2155090 | 5.5x90 mm  |
| 2155100 | 5.5x100 mm |
| 2155110 | 5.5x110 mm |
| 2155120 | 5.5x120 mm |
| 2155130 | 5.5x130 mm |
| 2155140 | 5.5x140 mm |
| 2155150 | 5.5x150 mm |
| 2155160 | 5.5x160 mm |
| 2155170 | 5.5x170 mm |
| 2155180 | 5.5x180 mm |
| 2155190 | 5.5x190 mm |
| 2155200 | 5.5x200 mm |

# Posterior Spinal Rods

## Features & Benefits

- Ti6Al4V-ELI raw material
- Easy screw fixation
- Flexibility for rod replacement
- Appropriate straight & curved rod in a variety of lengths
- Alternatively rods may be cut to length using the cutter
- The straight rod can be curved the desired lordosis using the Rod Bender
- 5.5 and 6.0 option



**Posterior Spinal  
Rods, 6.0**

### Spinal Rod 6.0

| Code    | Size       | Code    | Size       |
|---------|------------|---------|------------|
| 2200640 | 6.0x40 mm  | 2206270 | 6.0x270 mm |
| 2200650 | 6.0x50 mm  | 2206280 | 6.0x280 mm |
| 2200660 | 6.0x60 mm  | 2206290 | 6.0x290 mm |
| 2200670 | 6.0x70 mm  | 2206300 | 6.0x300 mm |
| 2200680 | 6.0x80 mm  | 2206310 | 6.0x310 mm |
| 2200690 | 6.0x90 mm  | 2206320 | 6.0x320 mm |
| 2206100 | 6.0x100 mm | 2206340 | 6.0x340 mm |
| 2206110 | 6.0x110 mm | 2206350 | 6.0x350 mm |
| 2206120 | 6.0x120 mm | 2206360 | 6.0x360 mm |
| 2206130 | 6.0x130 mm | 2206380 | 6.0x380 mm |
| 2206140 | 6.0x140 mm | 2206400 | 6.0x400 mm |
| 2206150 | 6.0x150 mm | 2206420 | 6.0x420 mm |
| 2206160 | 6.0x160 mm | 2206440 | 6.0x440 mm |
| 2206170 | 6.0x170 mm | 2206450 | 6.0x450 mm |
| 2206180 | 6.0x180 mm | 2206460 | 6.0x460 mm |
| 2206190 | 6.0x190 mm | 2206480 | 6.0x480 mm |
| 2206200 | 6.0x200 mm | 2206500 | 6.0x500 mm |
| 2206210 | 6.0x210 mm | 2206520 | 6.0x520 mm |
| 2206220 | 6.0x220 mm | 2206540 | 6.0x540 mm |
| 2206230 | 6.0x230 mm | 2206550 | 6.0x550 mm |
| 2206240 | 6.0x240 mm | 2206560 | 6.0x560 mm |
| 2206250 | 6.0x250 mm | 2206580 | 6.0x580 mm |
| 2206260 | 6.0x260 mm | 2206600 | 6.0x600 mm |

### 4U SPINAL HIGH FLEX ROD 6.0

| Code    | Size       |
|---------|------------|
| 2100635 | 6.0x35 mm  |
| 2100640 | 6.0x40 mm  |
| 2100645 | 6.0x45 mm  |
| 2100650 | 6.0x50 mm  |
| 2100655 | 6.0x55 mm  |
| 2100660 | 6.0x60 mm  |
| 2100665 | 6.0x65 mm  |
| 2100670 | 6.0x70 mm  |
| 2100680 | 6.0x80 mm  |
| 2100690 | 6.0x90 mm  |
| 2106100 | 6.0x100 mm |
| 2106110 | 6.0x110 mm |
| 2106120 | 6.0x120 mm |
| 2106130 | 6.0x130 mm |
| 2106140 | 6.0x140 mm |
| 2106150 | 6.0x150 mm |
| 2106160 | 6.0x160 mm |
| 2106170 | 6.0x170 mm |
| 2106180 | 6.0x180 mm |
| 2106190 | 6.0x190 mm |
| 2106200 | 6.0x200 mm |

# Thoracolumbar Connectors

## Multiaxial Transverse Connector

| Code    | Size  |
|---------|-------|
| 2600040 | 40 mm |
| 2600050 | 50 mm |
| 2600060 | 60 mm |
| 2600070 | 70 mm |
| 2600080 | 80 mm |



## Transverse Connector

| Code    | Size   |
|---------|--------|
| 2500040 | 40 mm  |
| 2500050 | 50 mm  |
| 2500060 | 60 mm  |
| 2500070 | 70 mm  |
| 2500080 | 80 mm  |
| 2500090 | 90 mm  |
| 2500100 | 100 mm |
| 2500120 | 120 mm |



## Lateral Connector

| Code    | Size  |
|---------|-------|
| 2700015 | 15 mm |
| 2700020 | 20 mm |
| 2700025 | 25 mm |
| 2700030 | 30 mm |



## Domino Connector

| Code    |
|---------|
| 2700040 |



## Axial Connector

| Code |
|------|
|------|

2700030



## Multiaxial Connector

| Code    | Size  |
|---------|-------|
| 2700050 | 50 mm |
| 2700060 | 60 mm |



# Thoracolumbar Hooks

## Pedicular Hooks Small

Code

Size

3000085  
3000087  
3000089




8x5  
8x7  
8x9

## Pedicular Offset Hooks Right

Code

Size

3100055  
3100057




5x5  
5x7

## Pedicular Offset Hooks Left

Code

Size

3200055  
3200057



5x5  
5x7

## Spinal Laminar Hooks

Code

Size

3300055  
3300057  
3300059  
3300076  
3300077  
3300079  
3300711




5x5  
5x7  
5x9  
7x6  
7x7  
7x9  
7x11

## Spinal Laminar Offset Hooks Left

Code

Size

3600055  
3600057




5x5  
5x7

## Spinal Laminar Offset Hooks Right

Code

Size

3500055  
3500057



5x5  
5x7

## Spinal Laminar Hooks Angled Left

Code

Size

3350079  
3350711



7x9  
7x11

## Spinal Laminar Hooks Angled Right

Code

Size

3450079  
3450711




7x9  
7x11


# Thoracolumbar Corpectomy Cage

## Features & Benefits

- Ti6Al4V-ELI raw material
- One-easy locking mechanism and teethed and angled superior and inferior surface minimize the displacement and coming out risk
- Rotatable superior & inferior endplate
- Efficient graft space
- Anatomical design
- Angled superior & inferior surface option



| Code    | D  | CL  | OL  |
|---------|----|-----|-----|
| 5702534 | 20 | 25  | 34  |
| 5703244 | 20 | 32  | 44  |
| 5704258 | 20 | 42  | 58  |
| 5705684 | 20 | 56  | 84  |
| 5802329 | 24 | 23  | 29  |
| 5802840 | 24 | 28  | 40  |
| 5803854 | 24 | 38  | 54  |
| 5805276 | 24 | 52  | 76  |
| 5807411 | 24 | 74  | 114 |
| 5903044 | 28 | 30  | 44  |
| 5904258 | 28 | 42  | 58  |
| 5905684 | 28 | 56  | 84  |
| 5908213 | 28 | 82  | 130 |
| 5901014 | 28 | 100 | 140 |



## Angled

|         |    |     |     |    |
|---------|----|-----|-----|----|
| 5752534 | 20 | 25  | 34  | 6° |
| 5753244 | 20 | 32  | 44  | 6° |
| 5754258 | 20 | 42  | 58  | 6° |
| 5755684 | 20 | 56  | 84  | 6° |
| 5852329 | 24 | 23  | 29  | 6° |
| 5852840 | 24 | 28  | 40  | 6° |
| 5853854 | 24 | 38  | 54  | 6° |
| 5855276 | 24 | 52  | 76  | 6° |
| 5857411 | 24 | 74  | 114 | 6° |
| 5953044 | 28 | 30  | 44  | 6° |
| 5954258 | 28 | 42  | 58  | 6° |
| 5955684 | 28 | 56  | 84  | 6° |
| 5958213 | 28 | 82  | 130 | 6° |
| 5951014 | 28 | 100 | 140 | 6° |

# Thoracolumbar Interspinous Device

## Features & Benefits


- Ti6Al4V-ELI raw material
- Easy usage
- Anatomical design

| Code    | Size |
|---------|------|
| 8800008 | 08   |
| 8800010 | 10   |
| 8800012 | 12   |
| 8800014 | 14   |
| 8800016 | 16   |



4U Spinal Interspinous  
U Device

# Expendable PLIF PEEK Cage




4U Spinal  
Expendable PLIF  
Cage

## Features & Benefits

- Poly Ether Ether Ketone (PEEK), Ti6Al4V-ELI raw material
- Toothed superior & inferior surface offer good and anatomical fixation
- Titanium Alloy mechanism helps to be seen in Imaging Machines
- Large graft area
- Compatible with MRI and CT


| Code    | L  | W  | H  |
|---------|----|----|----|
| 4602407 | 24 | 10 | 7  |
| 4602408 | 24 | 10 | 8  |
| 4602409 | 24 | 10 | 9  |
| 4602410 | 24 | 10 | 10 |
| 4602411 | 24 | 10 | 11 |
| 4602412 | 24 | 10 | 12 |
| 4602413 | 24 | 10 | 13 |



# Thoracolumbar PLIF PEEK Cage



## Features & Benefits

- Poly Ether Ether Ketone (PEEK) raw material
- Toothed superior & inferior surface offer good and anatomical fixation
- Tantalum markers helps to be seen in Imaging Machines
- Large graft area
- Compatible with MRI and CT




## 4U Spinal PLIF Cage

| Code    | L  | W  | H  | Code    | L  | W  | H  | Angled |
|---------|----|----|----|---------|----|----|----|--------|
| 4502206 | 22 | 10 | 6  | 4552206 | 22 | 10 | 6  | 5°     |
| 4502207 | 22 | 10 | 7  | 4552207 | 22 | 10 | 7  | 5°     |
| 4502208 | 22 | 10 | 8  | 4552208 | 22 | 10 | 8  | 5°     |
| 4502209 | 22 | 10 | 9  | 4552209 | 22 | 10 | 9  | 5°     |
| 4502210 | 22 | 10 | 10 | 4552210 | 22 | 10 | 10 | 5°     |
| 4502211 | 22 | 10 | 11 | 4552211 | 22 | 10 | 11 | 5°     |
| 4502212 | 22 | 10 | 12 | 4552212 | 22 | 10 | 12 | 5°     |
| 4502406 | 24 | 10 | 6  | 4552214 | 22 | 10 | 14 | 5°     |
| 4502407 | 24 | 10 | 7  | 4552216 | 22 | 10 | 16 | 5°     |
| 4502408 | 24 | 10 | 8  | 4552406 | 24 | 10 | 6  | 5°     |
| 4502409 | 24 | 10 | 9  | 4552407 | 24 | 10 | 7  | 5°     |
| 4502410 | 24 | 10 | 10 | 4552408 | 24 | 10 | 8  | 5°     |
| 4502411 | 24 | 10 | 11 | 4552409 | 24 | 10 | 9  | 5°     |
| 4502412 | 24 | 10 | 12 | 4552410 | 24 | 10 | 10 | 5°     |
|         |    |    |    | 4552411 | 24 | 10 | 11 | 5°     |
|         |    |    |    | 4552412 | 24 | 10 | 12 | 5°     |



# Thoracolumbar TLIF PEEK Cage



## Features & Benefits

- Poly Ether Ether Ketone (PEEK), Ti6Al4V-ELI raw material
- Toothed superior & inferior surface offer good and anatomical fixation
- Titanium Alloy mechanism helps to be seen in Imaging Machines
- 4° angle of lordosis
- Large graft area
- Compatible with MRI and CT

## 4U Spinal TLIF Cage

| Code    | L  | W  | H  | Code    | L  | W  | H  | Angled |
|---------|----|----|----|---------|----|----|----|--------|
| 4102407 | 24 | 10 | 7  | 4152407 | 24 | 10 | 7  | 4°     |
| 4102408 | 24 | 10 | 8  | 4152408 | 24 | 10 | 8  | 4°     |
| 4102409 | 24 | 10 | 9  | 4152409 | 24 | 10 | 9  | 4°     |
| 4102410 | 24 | 10 | 10 | 4152410 | 24 | 10 | 10 | 4°     |
| 4102411 | 24 | 10 | 11 | 4152411 | 24 | 10 | 11 | 4°     |
| 4102412 | 24 | 10 | 12 | 4152412 | 24 | 10 | 12 | 4°     |
| 4102413 | 24 | 10 | 13 | 4152413 | 24 | 10 | 13 | 4°     |
| 4202807 | 28 | 10 | 7  | 4252807 | 28 | 10 | 7  | 4°     |
| 4202808 | 28 | 10 | 8  | 4252808 | 28 | 10 | 8  | 4°     |
| 4202809 | 28 | 10 | 9  | 4252809 | 28 | 10 | 9  | 4°     |
| 4202810 | 28 | 10 | 10 | 4252810 | 28 | 10 | 10 | 4°     |
| 4202811 | 28 | 10 | 11 | 4252811 | 28 | 10 | 11 | 4°     |
| 4202812 | 28 | 10 | 12 | 4252812 | 28 | 10 | 12 | 4°     |
| 4202813 | 28 | 10 | 13 | 4252813 | 28 | 10 | 13 | 4°     |
| 4303207 | 32 | 10 | 7  | 4353207 | 32 | 10 | 7  | 4°     |
| 4303208 | 32 | 10 | 8  | 4353208 | 32 | 10 | 8  | 4°     |
| 4303209 | 32 | 10 | 9  | 4353209 | 32 | 10 | 9  | 4°     |
| 4303210 | 32 | 10 | 10 | 4353210 | 32 | 10 | 10 | 4°     |
| 4303211 | 32 | 10 | 11 | 4353211 | 32 | 10 | 11 | 4°     |
| 4303212 | 32 | 10 | 12 | 4353212 | 32 | 10 | 12 | 4°     |
| 4303213 | 32 | 10 | 13 | 4353213 | 32 | 10 | 13 | 4°     |


# Spinal Cervical Posterior Polyaxial Screw Set




# Spinal Posterior Polyaxial Screw Set




# Spinal MIS Set




# Spinal Cervical Plate Cage Set




# Spinal Cervical PEEK Cage Set





# Spinal PLIF & TLIF PEEK Cage Set



# Spinal Corpectomy Cage Set



# Spinal Interspinous U-Device Set




MEDICAL  
4U®

Turkey



[www.4umedical.com.tr](http://www.4umedical.com.tr)

# Keeps You Moving



**4U Medical**

-  Turgut Ozal Mah. 2118.Cd.No:6/A  
Yenimahalle/Ankara-TURKIYE
-  +90 312 566 06 55
-  +90 532 677 42 32
-  [info@4umedical.com.tr](mailto:info@4umedical.com.tr)



[www.4umedical.com.tr](http://www.4umedical.com.tr)

